

Anti-Acot2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P49753

Item Number: CM0014622



Introduction

Rabbit polyclonal antibody against ACOT2, validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat. Suitable for mitochondrial fatty acid oxidation research.

[Learn More](#)

Parameter	Value
Product Name	ACOT2 Rabbit pAb
Remarks/Alias	MTE1; PTE2; CTE1A; PTE2A; CTE-IA; ZAP128; ACOT2
Species	Human
GeneID (Human)	10965
GeneID	10965
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 222-483 of human ACOT2 (NP_006812.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49753
Protein Weight	53kDa
Shipping	Ice bag