

Anti-Acot8 Rabbit Polyclonal Antibody For Wb, Elisa - O14734

Item Number: CM0017374

Introduction

ACOT8 rabbit polyclonal antibody validated for Western blot (WB) and ELISA in human samples. Targets acyl-CoA thioesterase 8 (ACOT8), a peroxisomal enzyme involved in fatty acid metabolism, bile acid processing, and host-virus interactions.

[Learn More](#)

Parameter	Value
Product Name	ACOT8 Rabbit pAb
Remarks/Alias	hTE; NAP1; PTE1; PTE2; PTE-1; PTE-2; HNAACTE; hACTE-III; ACOT8
Species	Human
GeneID (Human)	10005
GeneID	10005
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ACOT8 (NP_005460.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14734
Protein Weight	36kDa
Shipping	Ice bag