

Anti-Acox1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15067

Item Number: CM0007482

Introduction

High-quality ACOX1 Rabbit mAb for WB, IHC-P, IF/ICC, ELISA. Specific for human, mouse, rat peroxisomal acyl-CoA oxidase 1, MW 74kDa. Ideal for metabolic and peroxisomal disorder research.

[Learn More](#)

Parameter	Value
Product Name	ACOX1 Rabbit mAb
Remarks/Alias	AOX, ACOX, SCOX, MITCH, PALMCOX, ACOX1
Species	Human
GeneID	51
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human ACOX1 (NP_004026.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15067
Protein Weight	74kDa
Shipping	Ice bag