

Anti-Acox1 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - Q15067

Item Number: CM0006573

Introduction

Rabbit monoclonal ACOX1 antibody validated for WB, IF/ICC, IF-P, IHC-P, ELISA with cross reactivity in human, mouse, rat. Ideal for peroxisomal fatty acid oxidation research.

[Learn More](#)

Parameter	Value
Product Name	ACOX1 Rabbit mAb
Remarks/Alias	AOX; ACOX; SCOX; MITCH; PALMCOX; ACOX1
Species	Human
Gene ID (Human)	51
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human ACOX1 (NP_004026.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15067
Protein Weight	74kDa
Shipping	Ice bag