

Anti-Acox2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99424

Item Number: CM0018399

Introduction

High-quality rabbit polyclonal antibody against human ACOX2, validated for WB, IHC-P, ELISA. Suitable for studies on peroxisomal beta-oxidation of branched-chain fatty acids and bile acid synthesis. Cross-reacts with human and mouse.

[Learn More](#)

Parameter	Value
Product Name	ACOX2 Rabbit pAb
Remarks/Alias	BCOX; BRCOX; CBAS6; THCCox; BRCACOX; ACOX2
Species	Human
Gene ID (Human)	8309
Gene ID	8309
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human ACOX2 (NP_003491.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99424
Protein Weight	77kDa
Shipping	Ice bag