

Anti-Acox3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O15254

Item Number: CM0020660

Introduction

Rabbit polyclonal antibody specifically recognizing peroxisomal acyl-coenzyme A oxidase 3 (ACOX3). Validated for WB, IF/ICC, ELISA applications in human, mouse, and rat samples. Ideal for branched-chain fatty acid metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ACOX3 Rabbit pAb
Remarks/Alias	ACOX3
Species	Human
Gene ID (Human)	8310
Gene ID	8310
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-100 of human ACOX3 (NP_001095137)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15254
Protein Weight	78kDa
Shipping	Ice bag