

Anti-Acp1 (Lmw-Ptp) Rabbit Polyclonal Antibody For Wb And Elisa - P24666

Item Number: CM0019546



Introduction

Rabbit polyclonal antibody against human ACP1 (LMW-PTP, HAAP). Validated for WB and ELISA, reacts with human and mouse. Suitable for studying phosphotyrosine protein phosphatase activity and related signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	ACP1 Rabbit pAb
Remarks/Alias	HAAP; LMWPTP; LMW-PTP; ACP1
Species	Human
GeneID (Human)	52
GeneID	52
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 15-144 of human ACP1 (NP_004291.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P24666
Protein Weight	18kDa
Shipping	Ice bag