


Anti-Acp5 Polyclonal Antibody For Wb, If, Ihc, Elisa - P13686

Item Number: CM0024112



Introduction

Rabbit polyclonal antibody targeting human tartrate-resistant acid phosphatase type 5 (ACP5/TRAP). Validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Cross-reacts with mouse and rat. Ideal for bone metabolism, osteoclast research, and leukemia biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	ACP5 Rabbit pAb
Remarks/Alias	HPAP; TRAP; TRAcP; TRACP5a; TRACP5b; TrATPase; ACP5
Species	Human
Gene ID (Human)	54
Gene ID	54
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 76-325 of human ACP5 (NP_001602.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13686
Protein Weight	37kDa
Shipping	Ice bag