

Anti-Acrv1 Polyclonal Antibody For Wb, Elisa - P26436

Item Number: CM0026288

Introduction

Rabbit polyclonal antibody targeting human ACRV1 (SP-10), an acrosomal protein essential for sperm function. Validated for Western blot and ELISA, cross-reacts with rat. Ideal for reproductive biology and infertility diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	ACRV1 Rabbit pAb
Remarks/Alias	SP-10, SPACA2, D11S4365, ACRV1
Species	Human
Gene ID (Human)	56
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-265 of human ACRV1 (NP_001603.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P26436
Protein Weight	28kDa
Shipping	Ice bag