

Anti-Acsf3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q4G176

Item Number: CM0019415

Introduction

Rabbit polyclonal antibody against human ACSF3 (malonate--CoA ligase), validated for WB, IHC-P, ELISA. Cross-reacts with mouse. Suitable for intramitochondrial fatty acid synthesis research.

[Learn More](#)

Parameter	Value
Product Name	ACSF3 Rabbit pAb
Remarks/Alias	ACSF3; acyl-CoA synthetase family member 3
Species	Human
Gene ID (Human)	197322
Gene ID	197322
Immunogen	Recombinant fusion protein containing amino acids 370-576 of human ACSF3 (NP_777577.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q4G176
Protein Weight	64kDa
Shipping	Ice bag