

Anti-Acs11 Polyclonal Antibody For Wb, Ip, If/Icc - P33121

Item Number: CM0020257

Introduction

High-quality rabbit polyclonal antibody against ACSL1, validated for WB, IP, IF/ICC in human, mouse, and rat samples. Ideal for fatty acid metabolism and lipid synthesis research.

[Learn More](#)

Parameter	Value
Product Name	ACSL1 Rabbit pAb
Remarks/Alias	ACS1; LACS; FACL1; FACL2; LACS1; LACS2; ACSL1
Species	Human
GeneID	2180
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 46-145 of human ACSL1 (NP_001986.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33121
Protein Weight	78 kDa
Shipping	Ice bag