

Anti-Acs11 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P33121

Item Number: CM0012295

Introduction

High-specificity rabbit monoclonal antibody targeting human ACSL1. Validated for WB, IF-P, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for studying long-chain fatty acid metabolism and lipid synthesis.

[Learn More](#)

Parameter	Value
Product Name	ACSL1 Rabbit mAb
Remarks/Alias	ACS1, LACS, FACL1, FACL2, LACS1, LACS2
Species	Human
Gene ID	2180
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 598-698 of human ACSL1 (NP_001986.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33121
Protein Weight	78kDa
Shipping	Ice bag