

Anti-Acs11 Rabbit Polyclonal Antibody For Wb, Elisa - P33121

Item Number: CM0021610

Introduction

High-quality rabbit polyclonal antibody targeting human ACSL1 (long-chain-fatty-acid--CoA ligase 1). Validated for Western blotting and ELISA, with cross-reactivity to mouse. Ideal for research into fatty acid activation, lipid metabolism, and metabolic disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	ACSL1 Rabbit pAb
Remarks/Alias	ACSL1; ACS1; FACL1; FACL2; LACS; LACS1; LACS2; long-chain-fatty-acid--CoA ligase 1
Species	Human
GeneID (Human)	2180
GeneID	2180
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P33121
Protein Weight	76kDa/77kDa
Shipping	Ice bag