

Anti-Acsl1 Rabbit Polymab® Monoclonal Antibody For Wb, Ihc, Elisa - P33121

Item Number: CM0007553

Introduction

ACSL1 Rabbit PolymAb® is a monoclonal antibody validated for WB, IF/ICC, IHC-P, and ELISA. Detects human, mouse, rat ACSL1, a long-chain fatty acid-CoA ligase involved in lipid metabolism. Ideal for metabolic research.

[Learn More](#)

Parameter	Value
Product Name	ACSL1 Rabbit PolymAb®
Remarks/Alias	ACS1; LACS; FACL1; FACL2; LACS1; LACS2
Species	Human
GeneID (Human)	2180
GeneID	2180
Immunogen	Synthetic peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ACSL1 (NP_001986.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33121
Protein Weight	78kDa
Shipping	Ice bag