

Anti-Acsl3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O95573

Item Number: CM0003968

Introduction

Rabbit monoclonal antibody against ACSL3, validated for WB, IF/ICC, ELISA. Reacts with human, mouse, and rat. Detects an 80kDa protein. Targets fatty acid CoA ligase for lipid metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ACSL3 Rabbit mAb
Remarks/Alias	ACS3; FACL3; LACS3; LACS 3; PRO2194; ACSL3
Species	Human
GeneID (Human)	2181
GeneID	2181
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 600-700 of human ACSL3 (NP_004448.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95573
Protein Weight	80kDa
Shipping	Ice bag