

Anti-Acsl4 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - O60488

Item Number: CM0002977

Introduction

Anti-ACSL4 Rabbit Monoclonal Antibody validated for WB, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat ACSL4 (O60488), a key enzyme in fatty acid metabolism and ferroptosis.

[Learn More](#)

Parameter	Value
Product Name	ACSL4 Rabbit mAb
Remarks/Alias	ACS4; FACL4; LACS4; MRX63; MRX68; XLID63; ACSL4
Species	Human
Gene ID (Human)	2182
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-321 of human ACSL4 (NP_075266.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60488
Protein Weight	79kDa
Shipping	Ice bag