

Anti-Acsl4 Rabbit Polyclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - O60488

Item Number: CM0026993

Introduction

Rabbit polyclonal antibody against human ACSL4, validated for WB, IP, IF/ICC, IHC-P, ELISA. Cross-reacts with Human, Mouse, Rat. Targets ACSL4, a key enzyme in long-chain fatty acid activation and ferroptosis.

[Learn More](#)

Parameter	Value
Product Name	ACSL4 Rabbit pAb
Remarks/Alias	ACS4; FACL4; LACS4; MRX63; MRX68; XLID63; ACSL4
Species	Human
GeneID (Human)	2182
GeneID	2182
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-321 of human ACSL4 (NP_075266.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60488
Protein Weight	79 kDa
Shipping	Ice bag