

Anti-AcsL5 Polyclonal Antibody For Wb And Elisa - Q9ULc5

Item Number: CM0017116

Introduction

Rabbit polyclonal antibody targeting human ACSL5, validated for Western blot and ELISA. Detects human, mouse, and rat ACSL5. Ideal for lipid metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ACSL5 Rabbit pAb
Remarks/Alias	ACS2; ACS5; FACL5; ACSL5
Species	Human
GeneID (Human)	51703
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 500-739 of human ACSL5 (NP_057318.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULC5
Protein Weight	76kDa
Shipping	Ice bag