

Anti-Acsl6 Polyclonal Antibody For Wb, Elisa - Q9Uku0

Item Number: CM0027059

Introduction

Rabbit polyclonal antibody against human ACSL6, validated for Western blot (WB) and ELISA. Detects endogenously expressed protein in human, mouse, and rat samples. Suitable for studies of long-chain fatty acid metabolism and brain lipid synthesis.

[Learn More](#)

Parameter	Value
Product Name	ACSL6 Rabbit pAb
Remarks/Alias	ACS2; FACL6; LACS2; LACS5; LACS 6; ACSL6
Species	Human
GeneID (Human)	23305
GeneID	23305
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 383-622 of human ACSL6 (NP_001192180.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKU0
Protein Weight	78kDa
Shipping	Ice bag