


Anti-Acsm2A Polyclonal Antibody For Wb, Elisa - Q08Ah3

Item Number: CM0014125



Introduction

Rabbit polyclonal antibody targeting human ACSM2A, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial fatty acid metabolism and xenobiotic activation research.

[Learn More](#)

Parameter	Value
Product Name	ACSM2A Rabbit pAb
Remarks/Alias	ACSM2; A-923A4.1; ACSM2A
Species	Human
GeneID (human)	123876
GeneID	123876
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 418-577 of human ACSM2A (NP_001295101.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q08AH3
Protein Weight	64kDa
Shipping	Ice bag