

## Anti-Acsm3 Rabbit Polyclonal Antibody For Wb, Elisa - Q53Fz2

Item Number: CM0019788

### Introduction

Rabbit polyclonal antibody targeting human ACSM3, validated for WB and ELISA. Cross-reacts with mouse. Useful for fatty acid metabolism research and mitochondrial studies.

[Learn More](#)

| Parameter        | Value                 |
|------------------|-----------------------|
| Product Name     | ACSM3 Rabbit pAb      |
| Remarks/Alias    | SA; SAH; ACSM3        |
| Species          | Human                 |
| Gene ID          | 6296                  |
| Immunogen        | synthetic peptide     |
| Source           | Rabbit                |
| Category         | Polyclonal Antibodies |
| Application      | WB, ELISA             |
| Cross Reactivity | Mouse                 |
| SWISS            | Q53FZ2                |
| Protein Weight   | 66kDa                 |
| Shipping         | Ice bag               |