


Anti-Acss1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nub1

Item Number: CM0013785



Introduction

Rabbit anti-ACSS1 polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat ACSS1. For acetyl-CoA metabolism, ketogenesis, and mitochondrial research. Ideal for studying acetate utilization, short-chain fatty acid metabolism, and thermogenesis.

[Learn More](#)

Parameter	Specification
Product Name	ACSS1 Rabbit pAb
Remarks/Alias	ACAS2L; ACECS1; AceCS2L; ACSS1
Species	Human
GeneID (Human)	84532
GeneID	84532
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-689 of human ACSS1 (NP_115890.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NUB1
Protein Weight	75kDa
Shipping	Ice bag