

## Anti-Acss2 Polyclonal Antibody For Wb, Elisa - Q9Nr19

Item Number: CM0016114

### Introduction

ACSS2 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and rat, targeting acetyl-CoA synthetase 2 (~79 kDa). Ideal for metabolic and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ACSS2 Rabbit pAb                                                                                                                        |
| Remarks/Alias    | ACS, ACSA, ACAS2, ACECS, AceCS1, dj1161H23.1, ACSS2                                                                                     |
| Species          | Human                                                                                                                                   |
| GeneID (Human)   | 55902                                                                                                                                   |
| GeneID           | 55902                                                                                                                                   |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 462-701 of human ACSS2 (NP_061147.1). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, ELISA                                                                                                                               |
| Cross Reactivity | Human, Rat                                                                                                                              |
| SWISS            | Q9NR19                                                                                                                                  |
| Protein Weight   | 79kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                 |