

Anti-Acss2 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Nr19

Item Number: CM0002631

Introduction

Validated rabbit monoclonal antibody against human ACSS2 (acetyl-CoA synthetase 2) for WB and ELISA applications. Detects mouse and rat ACSS2. Ideal for studying metabolic regulation, autophagy, and cancer epigenetics.

[Learn More](#)

Parameter	Value
Product Name	ACSS2 Rabbit mAb
Remarks/Alias	ACS; ACSA; ACAS2; ACECS; AceCS1; dJ1161H23.1; ACSS2
Species	Human
Gene ID (human)	55902
Gene ID	55902
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human ACSS2 (Q9NR19).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NR19
Protein Weight	79kDa
Shipping	Ice bag