

Anti-Acss2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nr19

Item Number: CM0026350

Introduction

High-quality rabbit polyclonal antibody against human ACSS2 (acetyl-CoA synthetase 2), validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ACSS2 Rabbit pAb
Remarks/Alias	ACS; ACSA; ACAS2; ACECS; AceCS1; dj1161H23.1; ACSS2
Species	Human
GeneID (Human)	55902
GeneID	55902
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 462-701 of human ACSS2 (NP_061147.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NR19
Protein Weight	79kDa
Shipping	Ice bag