


Anti-Activator Protein 2 (Ap-2/Tfap2A) Polyclonal Antibody For Wb, Elisa - P05549

Item Number: CM0016861



Introduction

Rabbit polyclonal antibody targeting human transcription factor AP-2-alpha (TFAP2A). Validated for WB and ELISA. Ideal for developmental biology and transcription regulation research.

[Learn More](#)

Parameter	Specification
Product Name	Activator protein 2 (AP-2/TFAP2A) Rabbit pAb
Remarks/Alias	AP-2; BOFS; AP2TF; TFAP2; AP-2alpha; Activator protein 2 (AP-2/TFAP2A)
Species	Human
GeneID (human)	7020
GeneID	7020
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-200 of human Activator protein 2 (AP-2/TFAP2A) (NP_003211.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P05549
Protein Weight	48kDa
Shipping	Ice bag