


Anti-Active Caspase-3 Monoclonal Antibody For Wb, Elisa - P42574

Item Number: CM0002467



Introduction

Rabbit monoclonal antibody targeting active Caspase-3 (CPP32, SCA-1). Validated for WB and ELISA in human samples. Ideal for apoptosis research. Reacts with human. Molecular weight 32 kDa.

[Learn More](#)

Parameter	Value
Product Name	active Caspase-3 Rabbit mAb
Remarks/Alias	CPP32; SCA-1; CPP32B; Active Caspase-3
Species	Human
GeneID (human)	836
GeneID	836
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CASP3 (P42574).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P42574
Protein Weight	32kDa
Shipping	Ice bag