


Anti-Actl6A Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O96019

Item Number: CM0012231



Introduction

Rabbit monoclonal antibody against human ACTL6A (BAF53A), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for chromatin remodeling, neurogenesis, and DNA repair research.

[Learn More](#)

Parameter	Specification
Product Name	ACTL6A Rabbit mAb
Remarks/Alias	Arp4; ACTL6; BAF53A; INO80K; SMARCN1; ARPN-BETA; ACTL6A
Species	Human
GeneID (Human)	86
GeneID	86
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 201-290 of human ACTL6A (O96019)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O96019
Protein Weight	47kDa
Shipping	Ice bag