


Anti-Actl6A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O96019

Item Number: CM0026610



Introduction

High-quality rabbit polyclonal antibody targeting ACTL6A (BAF53A), validated for WB, IHC-P, ELISA with human, mouse, rat cross-reactivity. Ideal for chromatin remodeling and neurogenesis research.

[Learn More](#)

Parameter	Value
Product Name	ACTL6A Rabbit pAb
Remarks/Alias	Arp4; ACTL6; BAF53A; INO80K; SMARCN1; ARPN-BETA; ACTL6A
Species	Human
Gene ID (Human)	86
Gene ID	86
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-120 of human ACTL6A (NP_004292.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O96019
Protein Weight	47kDa
Shipping	Ice bag