

Anti-Actl6B Polyclonal Antibody (Baf53B) For Wb/Ihc/If/Elisa - O94805

Item Number: CM0025789

Introduction

High-quality rabbit polyclonal antibody against human ACTL6B (BAF53B), validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for chromatin remodeling and neurodevelopment research.

[Learn More](#)

Parameter	Value
Product Name	ACTL6B Rabbit pAb
Remarks/Alias	ACTL6; DEE76; BAF53B; EIEE76; IDDSSAD; SMARCN2; arpNalpha; ACTL6B
Species	Human
Gene ID (Human)	51412
Gene ID	51412
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ACTL6B (NP_057272.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94805
Protein Weight	47kDa
Shipping	Ice bag