


Anti-Actn1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P12814

Item Number: CM0017109



Introduction

High-quality rabbit polyclonal antibody targeting human alpha-actinin-1 (ACTN1). Validated for WB, IHC-P, and ELISA; cross-reacts with human and mouse. Ideal for cytoskeleton, cell adhesion, and T-cell activation studies.

[Learn More](#)

Parameter	Value
Product Name	ACTN1 Rabbit pAb
Remarks/Alias	BDPLT15; ACTN1
Species	Human
GeneID (human)	87
GeneID	87
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human ACTN1 (NP_001123477.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P12814
Protein Weight	103kDa
Shipping	Ice bag