


Anti-Actn2 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P35609

Item Number: CM0005516



Introduction

Anti-ACTN2 rabbit monoclonal antibody validated for Western blot, immunofluorescence, and ELISA in human, mouse, and rat samples. Ideal for studying sarcomere structure and muscle diseases.

[Learn More](#)

Parameter	Value
Product Name	ACTN2 Rabbit mAb
Remarks/Alias	MPD6, CMH23, CMYP8, CMD1AA, MYOC0Z, ACTN2
Species	Human
Gene ID (Human)	88
Gene ID	88
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ACTN2 (P35609)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35609
Protein Weight	104kDa
Shipping	Ice bag