


Anti-Actn2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P35609

Item Number: CM0028076



Introduction

High-quality rabbit polyclonal antibody against human ACTN2 (alpha-actinin-2). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, and rat samples. Ideal for muscle biology and cardiomyopathy research.

[Learn More](#)

Parameter	Value
Product Name	ACTN2 Rabbit pAb
Remarks/Alias	MPD6; CMH23; CMYP8; CMD1AA; MYOCOZ; ACTN2
Species	Human
GeneID (human)	88
GeneID	88
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ACTN2 (NP_001094.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35609
Protein Weight	104 kDa
Shipping	Ice bag