

Anti-Actn3 Monoclonal Antibody For Wb, If-P, Elisa - Q08043

Item Number: CM0003950



Introduction

Recombinant rabbit monoclonal antibody targeting alpha-actinin-3 (ACTN3), suitable for WB, IF-P, ELISA. Cross-reacts with mouse and rat; detects a 103 kDa protein. Useful for fast skeletal muscle fiber analysis.

[Learn More](#)

Parameter	Value
Product Name	ACTN3 Rabbit mAb
Remarks/Alias	ACTN3D; ACTN3
Species	Human
GeneID (Human)	89
GeneID	89
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-800 of human ACTN3 (Q08043)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q08043
Protein Weight	103kDa
Shipping	Ice bag