

Anti-Actn3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q08043

Item Number: CM0018607



Introduction

High-quality rabbit polyclonal anti-ACTN3 antibody for Western blot, IHC-P, and ELISA. Detects alpha-actinin-3 in human, mouse, rat; explores fast skeletal muscle fiber actin cross-linking.

[Learn More](#)

Parameter	Value
Product Name	ACTN3 Rabbit pAb
Remarks/Alias	ACTN3D; ACTN3
Species	Human
GeneID (Human)	89
GeneID	89
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 580-630 of human ACTN3 (NP_001095.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08043
Protein Weight	103kDa
Shipping	Ice bag