

Anti-Actr10 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nz32

Item Number: CM0017020



Introduction

Rabbit polyclonal antibody targeting human ACTR10 (Arp11) for WB and ELISA applications. Cross-reacts with mouse and rat. Suitable for dynein/dynactin studies.

[Learn More](#)

Parameter	Value
Product Name	ACTR10 Rabbit pAb
Remarks/Alias	Arp10; Arp11; ACTR11; HARP11; ACTR10
Species	Human
GeneID (Human)	55860
GeneID	55860
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 218-417 of human ACTR10 (NP_060947.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZ32
Protein Weight	46kDa
Shipping	Ice bag