


Anti-Actr1A Rabbit Polyclonal Antibody For Wb, Elisa - P61163

Item Number: CM0030579



Introduction

Rabbit anti-ACTR1A polyclonal antibody optimized for WB and ELISA. Validated across Human, Mouse, Rat. Targets alpha-centractin, a dynactin complex component that activates dynein motor transport along microtubules. Ideal for cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	ACTR1A Rabbit pAb
Remarks/Alias	ARP1; Arp1A; CTRN1; ACTR1A
Species	Human
GeneID (Human)	10121
GeneID	10121
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 187-376 of human ACTR1A (NP_005727.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61163
Protein Weight	43kDa
Shipping	Ice bag