

Anti-Actr1B Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P42025

Item Number: CM0007470

Introduction

Rabbit monoclonal antibody targeting human ACTR1B (beta-centractin), validated for WB, IHC-P, IP, ELISA. Cross-reacts with rat. Ideal for centrosome and cytoskeleton research, and microtubule-based motility studies.

[Learn More](#)

Parameter	Value
Product Name	ACTR1B Rabbit mAb
Remarks/Alias	PC3; ARP1B; CTRN2; ACTR1B
Species	Human
GenelD (Human)	10120
GenelD	10120
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	P42025
Protein Weight	42kDa
Shipping	Ice bag