

Anti-Actr3B Polyclonal Antibody For Wb And Elisa - Q9P1U1

Item Number: CM0027199

Introduction

ACTR3B rabbit polyclonal antibody for WB and ELISA applications in human and mouse samples. Target functions in actin cytoskeleton organization and may impact tumor metastasis. Research use only.

[Learn More](#)

Parameter	Value
Product Name	ACTR3B Rabbit pAb
Remarks/Alias	ARP11; ARP3BETA; ACTR3B
Species	Human
GeneID (Human)	57180
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 189-418 of human ACTR3B (NP_065178.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P1U1
Protein Weight	48kDa
Shipping	Ice bag