

Anti-Actrt3 Polyclonal Antibody For Wb, Elisa - Q9Byd9

Item Number: CM0016980

Introduction

Rabbit polyclonal antibody against human ACTRT3, validated for WB and ELISA. Cross-reacts with mouse and rat, detects a 41 kDa band. For cytoskeleton and nuclear protein research.

[Learn More](#)

Parameter	Value
Product Name	ACTRT3 Rabbit pAb
Remarks/Alias	ARPM1, ARP-T3, ACTRT3
Species	Human
Gene ID (Human)	84517
Gene ID	84517
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human ACTRT3 (NP_115876.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYD9
Protein Weight	41kDa
Shipping	Ice bag