


Anti-Acvr1 Rabbit Polyclonal Antibody (Pab) For Wb, Elisa - Q04771

Item Number: CM0026869



Introduction

Rabbit polyclonal antibody against human ACVR1 (activin receptor type-1, ALK2) for western blot (WB) and ELISA. Immunogen: recombinant protein (aa 20-125). Detects human and mouse ACVR1. Ideal for BMP signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	ACVR1 Rabbit pAb
Remarks/Alias	FOP; ALK2; SKR1; TSRI; ACTRI; ACVR1A; ACVRLK2; ACVR1
Species	Human
GeneID (Human)	90
GeneID	90
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-125 of human ACVR1 (NP_001096.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q04771
Protein Weight	57kDa
Shipping	Ice bag