

Anti-Acvr1B Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P36896

Item Number: CM0006892



Introduction

High-quality rabbit monoclonal antibody against human ACVR1B (Activin receptor type-1B, ALK4). Validated for WB, IHC-P, IF/ICC, and ELISA applications. Shows cross-reactivity with mouse and rat. Suitable for activin signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	ACVR1B Rabbit mAb
Remarks/Alias	ALK4; SKR2; ACTRIB; ACVRLK4; ACVR1B
Species	Human
GeneID (human)	91
GeneID	91
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ACVR1B (P36896).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36896
Protein Weight	57kDa
Shipping	Ice bag