

Anti-Acvr1B Rabbit Polyclonal Antibody For Wb, Elisa - P36896

Item Number: CM0028931



Introduction

Rabbit polyclonal antibody against human ACVR1B (ALK4). Validated for WB and ELISA. Cross-reacts with human and mouse. Suitable for activin signaling and TGF- β superfamily research.

[Learn More](#)

Parameter	Value
Product Name	ACVR1B Rabbit pAb
Remarks/Alias	ALK4; SKR2; ACTR1B; ACVRLK4; ACVR1B
Species	Human
Gene ID (Human)	91
Gene ID	91
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-126 of human ACVR1B (NP_004293.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P36896
Protein Weight	57kDa
Shipping	Ice bag