

Anti-Acvr1C Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Ner5

Item Number: CM0014528

Introduction

Rabbit polyclonal antibody against ACVR1C (Activin receptor type-1C/ALK7). Tested in WB, IHC-P, ELISA with reactivity to human, mouse, rat. Suitable for research on TGF-beta signaling and pancreatic beta-cell function.

[Learn More](#)

Parameter	Value
Product Name	ACVR1C Rabbit pAb
Remarks/Alias	ALK7; ACVRLK7; ACVR1C
Species	Human
Gene ID (Human)	130399
Gene ID	130399
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-113 of human ACVR1C (NP_660302.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NER5
Protein Weight	55kDa
Shipping	Ice bag