

Anti-Acvr2A Rabbit Monoclonal Antibody For Wb And Elisa - P27037

Item Number: CM0012045



Introduction

Rabbit monoclonal antibody targeting human ACVR2A, activin receptor type-2A. Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for studying TGF-beta signaling pathways. - P27037.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ACVR2A Rabbit mAb
Remarks/Alias	ACVR2; ACTRII; 2A
Species	Human
GeneID (human)	92
GeneID	92
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ACVR2A (P27037)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27037
Protein Weight	58kDa
Shipping	Ice bag