

Anti-Acvr2A Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P27037

Item Number: CM0022995



Introduction

High-quality ACVR2A Rabbit pAb validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat ACVR2A (58 kDa). Suitable for investigating TGF- β /activin signaling and receptor biology.

[Learn More](#)

Parameter	Value
Product Name	ACVR2A Rabbit pAb
Remarks/Alias	ACVR2; ACTRII; ACVR2A
Species	Human
GeneID (Human)	92
GeneID	92
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27037
Protein Weight	58 kDa
Shipping	Ice bag