

Anti-Acvr2B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13705

Item Number: CM0026089



Introduction

High-quality rabbit polyclonal antibody targeting human ACVR2B, validated for WB, IF/ICC, and ELISA. Ideal for activin signaling research and diagnostic applications. Detects 58kDa protein.

[Learn More](#)

Specification	Detail
Product Name	ACVR2B Rabbit pAb
Remarks/Alias	HTX4; ACTRIIB; ActR-IIB; ACVR2B
Species	Human
Gene ID (Human)	93
Gene ID	93
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human ACVR2B (NP_001097.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q13705
Protein Weight	58kDa
Shipping	Ice bag