

Anti-Acvrl1 Monoclonal Antibody For Wb, If/Icc, Elisa - P37023

Item Number: CM0009825



Introduction

High-quality ACVRL1 (ALK1) rabbit monoclonal antibody for WB, IF/ICC, ELISA. Detects human, mouse, rat ACVRL1. Essential for studying blood vessel development and TGF-beta superfamily signaling.

[Learn More](#)

Parameter	Value
Product Name	ACVRL1 Rabbit mAb
Remarks/Alias	HHT; ALK1; HHT2; ORW2; SKR3; ALK-1; TSR-I; ACVRLK1; ACVRL1
Species	Human
GeneID (Human)	94
GeneID	94
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 404-503 of human ACVRL1 (P37023).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37023
Protein Weight	56kDa
Shipping	Ice bag