

Anti-Acvrl1/Alk1 Rabbit Polyclonal Antibody For Wb And Elisa - P37023

Item Number: CM0022703



Introduction

ACVRL1 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human, mouse, and rat. Recombinant immunogen targeting the intracellular domain of the TGF-beta receptor ALK1, a key angiogenesis regulator. 56 kDa.

[Learn More](#)

Parameter	Specification
Product Name	ACVRL1 Rabbit pAb
Remarks/Alias	HHT; ALK1; HHT2; ORW2; SKR3; ALK-1; TSR-I; ACVRLK1; ACVRL1
Species	Human
GeneID (Human)	94
GeneID	94
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 234-503 of human ACVRL1 (NP_001070869.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37023
Protein Weight	56kDa
Shipping	Ice bag