

Anti-Acy1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q03154

Item Number: CM0025343



Introduction

High-quality rabbit polyclonal antibody against ACY1, validated for WB, IF/ICC, ELISA.

Recognizes human, mouse, rat ACY1 (aminoacylase-1) involved in N-acylamino acid hydrolysis. Suitable for degradation pathway research.

[Learn More](#)

Parameter	Value
Product Name	ACY1 Rabbit pAb
Remarks/Alias	ACY-1; ACY1D; HEL-S-5; ACY1
Species	Human
GeneID (human)	95
GeneID	95
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-408 of human ACY1 (NP_001185824.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03154
Protein Weight	46kDa
Shipping	Ice bag